

## **ICP Multi-Element Standard - 250 mL 4 components mixture and 250 mL Hg, in 2% HNO<sub>3</sub>**

Art. ID LM24-CUS-70067  
Unit 2 x 250 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Conservation 8 months (4 month after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mixture A	Boron (B)	[7440-42-8]	1	mg/L		
Mixture A	Silicon (Si)	[7440-21-3]	1	mg/L		
Mixture A	Thallium (Tl)	[7440-28-0]	0,005	mg/L		
Mixture A	Uranium (U)	[7440-61-1]	0,01	mg/L		
Single solution	Mercury (Hg)	[7439-97-6]	0,001	mg/L		