

## **Proficiency Testing - Oxygen, dissolved in water**

Art. ID	ERA-212
Unit	500 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Applicable methods: SM 4500-O, EPA 360.1, EPA 360.1, ASTM D888-12e1 // The Dissolved Oxygen standard is provided as a 500 mL whole volume sample with no headspace in an amber glass bottle sealed in mylar packaging.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Oxygen, Dissolved		1,0 - 20,0	mg/L	SM 4500-O, EPA 360.1 , EPA 360. 1, ASTM D8 88-12e1	