

Custom Solution - 6 components in Methanol each 500 mg/L

Art. ID LM24-CUS-70068
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Aminophenylether (4,4'-Oxydianiline)	[101-80-4]	500	mg/L		
	2-Methoxy-5-methylaniline	[120-71-8]	500	mg/L		
	3,3'-Dimethylbenzidine (o-Tolidine)	[119-93-7]	500	mg/L		
	2,4,5-Trimethylaniline	[137-17-7]	500	mg/L		
	4-Chloro-2-methylaniline	[95-69-2]	500	mg/L		
	Naphthalene D8	[1146-65-2]	500	mg/L		