

Custom Solution - 11 components in 5% HNO₃

Art. ID LM24-CUS-70073
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Reference value	Zinc (Zn)	[7440-66-6]	100	mg/L		
Reference value	Copper (Cu)	[7440-50-8]	100	mg/L		
Reference value	Lithium (Li)	[7439-93-2]	100	mg/L		
Reference value	Mercury (Hg)	[7439-97-6]	100	mg/L		
Reference value	Arsenic (As)	[7440-38-2]	100	mg/L		
Reference value	Cadmium (Cd)	[7440-43-9]	100	mg/L		
Reference value	Lead (Pb)	[7439-92-1]	100	mg/L		
Reference value	Thallium (Tl)	[7440-28-0]	100	mg/L		
Reference value	Aluminium (Al)	[7429-90-5]	100	mg/L		
Reference value	Chromium (Cr)	[7440-47-3]	100	mg/L		
Reference value	Selenium (Se)	[7782-49-2]	100	mg/L		