

Scallop

Art. ID NIES-15
Unit 20 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|-------------|---------------------|-------|---|--------|
| Certified value | TBT [Sn(C ₄ H ₉) ₃ +] | | 0,404 ± 0,027 | mg/kg | GC-MS,LC-M S | |
| Certified value | TPT [Sn(C ₆ H ₅) ₃ +] | | 0,017 ± 0,0017 | mg/kg | GC-MS,LC-M S | |
| Certified value | Tin (Sn) | [7440-31-5] | 0,179 ± 0,021 | mg/kg | ICP-MS,ID- ICP-MS,INA A | |
| Reference value | Sodium (Na) | [7440-23-5] | 0,561 | % | ICP-AES,IN AA | |
| Reference value | Magnesium (Mg) | [7439-95-4] | 0,179 | % | ICP-AES,IN AA,XRF | |
| Reference value | Phosphorus (P) | [7723-14-0] | 0,956 | % | ICP-AES,XR F | |
| Reference value | Potassium (K) | [7440-09-7] | 1,65 | % | ICP-AES,IN AA,XRF | |
| Reference value | Iron (Fe) | [7439-89-6] | 10,8 | mg/kg | ICP-AES,IC P-MS | |
| Reference value | Copper (Cu) | [7440-50-8] | 0,889 | mg/kg | ICP-AES,IC P-MS | |
| Reference value | Zinc (Zn) | [7440-66-6] | 56,1 | mg/kg | ICP-AES,IC P-MS,NAA | |
| Reference value | Arsenic (As) | [7440-38-2] | 3,76 | mg/kg | AS,F-AAS,H R-ICP-MS,I CP-MS,INAA ,MIP-MS,XR F | |
| Reference value | Bromine (Br) | [7726-95-6] | 43,7 | mg/kg | INAA,XRF | |