

### **Custom Solution mixture- 3 components (deuterated) in Methanol**

Art. ID LM24-CUS-70079  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-butyl ester D4 (Bis-butyl phthalate D4)	[93952-11-5]	100	µg/mL		
Di-octyl phthalate-3,4,5,6-d4	Di-n-octyl-phthalate D4	[93952-13-7]	100	µg/mL		
	Phthalic acid, benzylbutyl ester (3,4,5,6-D4) (Benzylbutyl phthalate D4)	[93951-88-3]	100	µg/mL		