

2,4-D 1000 µg/mL in Acetonitrile CERTAN®

Art. ID AP-CE445-4.5A
Unit 4.5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | 2,4-D | [94-75-7] | 1000 | µg/mL | | |