

Multielement aqueous calibration solution in 2% HNO₃

Art. ID AP-CZ9082MN1
Unit 100 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Barium (Ba) | [7440-39-3] | 1000 ± 2 | mg/L | | |
| | Calcium (Ca) | [7440-70-2] | 1000 ± 2 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 1000 ± 2 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 1000 ± 2 | mg/L | | |