

### **Multielement aqueous calibration solution in 5% aqua regia**

Art. ID AP-CZ9085MCN1  
Unit 100 mL/glass  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	25 ± 0.05	mg/L		
	Beryllium (Be)	[7440-41-7]	10 ± 0.02	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10 ± 0.02	mg/L		
	Cadmium (Cd)	[7440-43-9]	2 ± 0.004	mg/L		
	Mercury (Hg)	[7439-97-6]	2 ± 0.004	mg/L		
	Cobalt (Co)	[7440-48-4]	50 ± 0.1	mg/L		
	Copper (Cu)	[7440-50-8]	50 ± 0.1	mg/L		
	Chromium (Cr)	[7440-47-3]	100 ± 0.2	mg/L		
	Nickel (Ni)	[7440-02-0]	100 ± 0.2	mg/L		
	Zinc (Zn)	[7440-66-6]	100 ± 0.2	mg/L		
	Lead (Pb)	[7439-92-1]	200 ± 0.4	mg/L		
	Vanadium (V)	[7440-62-2]	200 ± 0.4	mg/L		