

### **Multielement aqueous calibration solution in 5% HNO<sub>3</sub> + 0.2% HF**

Art. ID AP-CZ9094MFN1  
Unit 100 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100 ± 0.2	mg/L		
	Beryllium (Be)	[7440-41-7]	100 ± 0.2	mg/L		
	Calcium (Ca)	[7440-70-2]	100 ± 0.2	mg/L		
	Cadmium (Cd)	[7440-43-9]	100 ± 0.2	mg/L		
	Cobalt (Co)	[7440-48-4]	100 ± 0.2	mg/L		
	Chromium (Cr)	[7440-47-3]	100 ± 0.2	mg/L		
	Copper (Cu)	[7440-50-8]	100 ± 0.2	mg/L		
	Iron (Fe)	[7439-89-6]	100 ± 0.2	mg/L		
	Lithium (Li)	[7439-93-2]	100 ± 0.2	mg/L		
	Magnesium (Mg)	[7439-95-4]	100 ± 0.2	mg/L		
	Manganese (Mn)	[7439-96-5]	100 ± 0.2	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100 ± 0.2	mg/L		
	Nickel (Ni)	[7440-02-0]	100 ± 0.2	mg/L		
	Phosphorus (P)	[7723-14-0]	100 ± 0.2	mg/L		
	Lead (Pb)	[7439-92-1]	100 ± 0.2	mg/L		
	Antimony (Sb)	[7440-36-0]	100 ± 0.2	mg/L		
	Selenium (Se)	[7782-49-2]	100 ± 0.2	mg/L		
	Tin (Sn)	[7440-31-5]	100 ± 0.2	mg/L		
	Strontium (Sr)	[7440-24-6]	100 ± 0.2	mg/L		
	Titanium (Ti)	[7440-32-6]	100 ± 0.2	mg/L		
	Thallium (Tl)	[7440-28-0]	100 ± 0.2	mg/L		
	Vanadium (V)	[7440-62-2]	100 ± 0.2	mg/L		
	Zinc (Zn)	[7440-66-6]	100 ± 0.2	mg/L		