

Multielement aqueous calibration solution in 5% HNO₃

Art. ID AP-CZ9095MN1
Unit 100 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100 ± 0.2	mg/L		
	Beryllium (Be)	[7440-41-7]	100 ± 0.2	mg/L		
	Bismuth (Bi)	[7440-69-9]	100 ± 0.2	mg/L		
	Cadmium (Cd)	[7440-43-9]	100 ± 0.2	mg/L		
	Cobalt (Co)	[7440-48-4]	100 ± 0.2	mg/L		
	Copper (Cu)	[7440-50-8]	100 ± 0.2	mg/L		
	Iron (Fe)	[7439-89-6]	100 ± 0.2	mg/L		
	Lithium (Li)	[7439-93-2]	100 ± 0.2	mg/L		
	Manganese (Mn)	[7439-96-5]	100 ± 0.2	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100 ± 0.2	mg/L		
	Nickel (Ni)	[7440-02-0]	100 ± 0.2	mg/L		
	Lead (Pb)	[7439-92-1]	100 ± 0.2	mg/L		
	Silicon (Si)	[7440-21-3]	100 ± 0.2	mg/L		
	Strontium (Sr)	[7440-24-6]	100 ± 0.2	mg/L		
	Vanadium (V)	[7440-62-2]	100 ± 0.2	mg/L		
	Tungsten (W)	[7440-33-7]	100 ± 0.2	mg/L		
	Zinc (Zn)	[7440-66-6]	100 ± 0.2	mg/L		
	Zirconium (Zr)	[7440-67-7]	100 ± 0.2	mg/L		