

### **Multielement aqueous calibration solution in 5% HNO<sub>3</sub>**

Art. ID AP-CZ9096MN1  
Unit 100 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100 ± 0.2	mg/L		
	Antimony (Sb)	[7440-36-0]	100 ± 0.2	mg/L		
	Selenium (Se)	[7782-49-2]	100 ± 0.2	mg/L		
	Tin (Sn)	[7440-31-5]	100 ± 0.2	mg/L		