

Multielement aqueous calibration solution in 5% HNO₃

Art. ID AP-CZ9099MN1
Unit 100 mL/glass
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100 ± 0.2	mg/L		
	Mercury (Hg)	[7439-97-6]	100 ± 0.2	mg/L		