

## **Multielement aqueous calibration solution in Water**

Art. ID AP-CZ9103MH1  
Unit 100 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100 ± 0.2	mg/L		
	Calcium (Ca)	[7440-70-2]	100 ± 0.2	mg/L		
	Potassium (K)	[7440-09-7]	100 ± 0.2	mg/L		
	Lithium (Li)	[7439-93-2]	100 ± 0.2	mg/L		
	Magnesium (Mg)	[7439-95-4]	100 ± 0.2	mg/L		
	Sodium (Na)	[7440-23-5]	100 ± 0.2	mg/L		
	Ammonium (NH <sub>4</sub> <sup>+</sup> )	[14798-03-9]	100 ± 0.2	mg/L		
	Strontium (Sr)	[7440-24-6]	100 ± 0.2	mg/L		