

### **Aqueous calibration solution (uncertified) in 5% HNO<sub>3</sub>**

Art. ID AP-TUNE-1  
Unit 25 mL/HDPE  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10 ± 0.1	mg/L		
	Beryllium (Be)	[7440-41-7]	10 ± 0.1	mg/L		
	Bismuth (Bi)	[7440-69-9]	10 ± 0.1	mg/L		
	Cerium (Ce)	[7440-45-1]	10 ± 0.1	mg/L		
	Cobalt (Co)	[7440-48-4]	10 ± 0.1	mg/L		
	Indium (In)	[7440-74-6]	10 ± 0.1	mg/L		
	Lithium (Li)	[7439-93-2]	10 ± 0.1	mg/L		
	Nickel (Ni)	[7440-02-0]	10 ± 0.1	mg/L		
	Lead (Pb)	[7439-92-1]	10 ± 0.1	mg/L		
	Uranium (U)	[7440-61-1]	10 ± 0.1	mg/L		