

Aqueous calibration solution (uncertified) in 5% HNO₃

Art. ID AP-TUNE-2
Unit 25 mL/HDPE
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	10 ± 0.1	mg/L		
	Iron (Fe)	[7439-89-6]	10 ± 0.1	mg/L		
	Potassium (K)	[7440-09-7]	10 ± 0.1	mg/L		
	Lithium (Li)	[7439-93-2]	10 ± 0.1	mg/L		
	Sodium (Na)	[7440-23-5]	10 ± 0.1	mg/L		