

### **Aqueous calibration solution (uncertified) in 5% HNO<sub>3</sub>**

Art. ID AP-TUNE-3  
Unit 25 mL/HDPE  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10 ± 0.1	mg/L		
	Arsenic (As)	[7440-38-2]	10 ± 0.1	mg/L		
	Barium (Ba)	[7440-39-3]	10 ± 0.1	mg/L		
	Beryllium (Be)	[7440-41-7]	10 ± 0.1	mg/L		
	Calcium (Ca)	[7440-70-2]	10 ± 0.1	mg/L		
	Cadmium (Cd)	[7440-43-9]	10 ± 0.1	mg/L		
	Cobalt (Co)	[7440-48-4]	10 ± 0.1	mg/L		
	Chromium (Cr)	[7440-47-3]	10 ± 0.1	mg/L		
	Copper (Cu)	[7440-50-8]	10 ± 0.1	mg/L		
	Iron (Fe)	[7439-89-6]	10 ± 0.1	mg/L		
	Indium (In)	[7440-74-6]	10 ± 0.1	mg/L		
	Potassium (K)	[7440-09-7]	10 ± 0.1	mg/L		
	Lithium (Li)	[7439-93-2]	10 ± 0.1	mg/L		
	Magnesium (Mg)	[7439-95-4]	10 ± 0.1	mg/L		
	Manganese (Mn)	[7439-96-5]	10 ± 0.1	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10 ± 0.1	mg/L		
	Sodium (Na)	[7440-23-5]	10 ± 0.1	mg/L		
	Nickel (Ni)	[7440-02-0]	10 ± 0.1	mg/L		
	Lead (Pb)	[7439-92-1]	10 ± 0.1	mg/L		
	Thallium (Tl)	[7440-28-0]	10 ± 0.1	mg/L		
	Vanadium (V)	[7440-62-2]	10 ± 0.1	mg/L		
	Uranium (U)	[7440-61-1]	10 ± 0.1	mg/L		
	Zinc (Zn)	[7440-66-6]	10 ± 0.1	mg/L		