

Aqueous calibration solution (uncertified) in 5% HNO₃

Art. ID AP-TUNE-4
Unit 25 mL/HDPE
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	10 ± 0.1	mg/L		
	Indium (In)	[7440-74-6]	10 ± 0.1	mg/L		
	Lead (Pb)	[7439-92-1]	10 ± 0.1	mg/L		
	Selenium (Se)	[7782-49-2]	10 ± 0.1	mg/L		
	Vanadium (V)	[7440-62-2]	10 ± 0.1	mg/L		