

Aqueous calibration solution (uncertified) in 5% HNO₃

Art. ID AP-INT-MIX-1
Unit 25 mL/HDPE
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	10 ± 0.1	mg/L		
	Indium (In)	[7440-74-6]	10 ± 0.1	mg/L		
	Scandium (Sc)	[7440-20-2]	10 ± 0.1	mg/L		
	Terbium (Tb)	[7440-27-9]	10 ± 0.1	mg/L		
	Yttrium (Y)	[7440-65-5]	10 ± 0.1	mg/L		