

Density Standard 1.8366g/ml at 15°C (ASTM D1480)

Art. ID REA-DEN15225PY
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 15°C | Density | | 1,8366 | g/mL | | |