

### **Multielement aqueous calibration solution (uncertified) in 2% HNO<sub>3</sub>**

Art. ID AP-M1007  
Unit 500 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	200 ± 1	mg/L		
	Barium (Ba)	[7440-39-3]	200 ± 1	mg/L		
	Calcium (Ca)	[7440-70-2]	200 ± 1	mg/L		
	Chromium (Cr)	[7440-47-3]	200 ± 1	mg/L		
	Copper (Cu)	[7440-50-8]	200 ± 1	mg/L		
	Iron (Fe)	[7439-89-6]	200 ± 1	mg/L		
	Magnesium (Mg)	[7439-95-4]	200 ± 1	mg/L		
	Manganese (Mn)	[7439-96-5]	200 ± 1	mg/L		
	Nickel (Ni)	[7440-02-0]	200 ± 1	mg/L		
	Silicon (Si)	[7440-21-3]	200 ± 1	mg/L		
	Strontium (Sr)	[7440-24-6]	200 ± 1	mg/L		