

Multielement aqueous calibration solution (uncertified) in 2% HNO₃

Art. ID AP-M1012
Unit 500 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	1000 ± 5	mg/L		
	Calcium (Ca)	[7440-70-2]	500 ± 2.5	mg/L		
	Phosphorus (P)	[7723-14-0]	400 ± 2	mg/L		
	Sodium (Na)	[7440-23-5]	240 ± 1.2	mg/L		
	Magnesium (Mg)	[7439-95-4]	100 ± 0.5	mg/L		
	Iron (Fe)	[7439-89-6]	10 ± 0.05	mg/L		
	Zinc (Zn)	[7440-66-6]	6 ± 0.03	mg/L		
	Copper (Cu)	[7440-50-8]	1 ± 0.005	mg/L		
	Manganese (Mn)	[7439-96-5]	1 ± 0.005	mg/L		