

Multielement (ion) calibration solution in Water

Art. ID AP-CZ9103MH1IC
Unit 100 mL/HDPE
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100 ± 0,2	mg/L		
	Calcium (Ca)	[7440-70-2]	100 ± 0,2	mg/L		
	Potassium (K)	[7440-09-7]	100 ± 0,2	mg/L		
	Lithium (Li)	[7439-93-2]	100 ± 0,2	mg/L		
	Magnesium (Mg)	[7439-95-4]	100 ± 0,2	mg/L		
	Sodium (Na)	[7440-23-5]	100 ± 0,2	mg/L		
	Ammonium (NH ₄ ⁺)	[14798-03-9]	100 ± 0,2	mg/L		
	Strontium (Sr)	[7440-24-6]	100 ± 0,2	mg/L		