

### **PAH MIX-4 in Acetonitrile CERTAN**

Art. ID AP-CE004-1.5A  
Unit 1.5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzo(a)pyrene	[50-32-8]	100	µg/mL		
	Fluoranthene	[206-44-0]	200	µg/mL		