

PAH MIX-11 in Toluene CERTAN

Art. ID AP-CE011-4.5T
Unit 4,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Anthracene	[120-12-7]	20	µg/mL		
	Benzo(a)anthracene	[56-55-3]	25	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	15	µg/mL		
	Benzo(j)fluoranthene	[205-82-3]	10	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	10	µg/mL		
	7H-Benzo(c)fluorene	[205-12-9]	10	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	10	µg/mL		
	Benzo(a)pyrene	[50-32-8]	15	µg/mL		
	Chrysene	[218-01-9]	30	µg/mL		
	Cyclopenta(c,d)pyrene	[27208-37-3]	20	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	2,5	µg/mL		
	Dibenzo(a,e)pyrene	[192-65-4]	15	µg/mL		
	Dibenzo(a,h)pyrene	[189-64-0]	15	µg/mL		
	Dibenzo(a,i)pyrene	[189-55-9]	15	µg/mL		
	Dibenzo(a,l)pyrene	[191-30-0]	15	µg/mL		
	Fluoranthene	[206-44-0]	80	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	15	µg/mL		
	5-Methylchrysene	[3697-24-3]	10	µg/mL		
	Phenanthrene	[85-01-8]	200	µg/mL		
	Pyrene	[129-00-0]	80	µg/mL		