

Custom Solution - 3 components 2500 µg/mL each of in Methanol

Art. ID LM24-CUS-70568
Unit 3x1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Bromofluorobenzene	[460-00-4]	2500	µg/mL		
	1,2-Dichloroethane D4	[17060-07-0]	2500	µg/mL		
	Toluene D8	[2037-26-5]	2500	µg/mL		