

Custom Solution - 3 components 2000 mg/L each of in Methylene chloride (Dichloromethane)

Art. ID LM24-CUS-70574
Unit 3 x 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Conservation: 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,4-Dichlorobenzene D4	[3855-82-1]	2000	mg/L		
	Naphthalene D8	[1146-65-2]	2000	mg/L		
	Acenaphthene D10	[15067-26-2]	2000	mg/L		