

Custom Solution - Phthalic acid, bis-iso-octyl ester 100 mg/L in Methylene chloride (Dichloromethane)

Art. ID LM24-CUS-70575
Unit 3 x 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Conservation: 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-iso-octyl ester (Bis-iso-octyl phthalate)	[27554-26-3]	100	mg/L		