

**PCB 90 500 µg/mL in Iso-octane**

Art. ID LM24-OP-4010268  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

| Text/Information | Analyte/Parameter                            | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
|                  | 2,2',3,4',5-Pentachloro<br>biphenyl (PCB 90) | [68194-07-0] | 500                 | µg/mL |        |        |