

**PCB 109 500 µg/mL in Iso-octane**

Art. ID LM24-OP-4010272  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

| Text/Information | Analyte/Parameter                        | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
|                  | 2,3,3',4,6-Pentachlorobiphenyl (PCB 109) | [74472-35-8] | 500                 | µg/mL |        |        |