

Pesticide Residues (2341)-Method IV GC/MS/MS :Mixed reference stock solution: Part 3 - 25 components, 1000 µg/L in Acetonitrile LC-MS

Art. ID LM24-OP-4010412
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------|---------------|---------------------|-------|--------|--------|
| cis-Bifenthrin | Bifenthrin | [82657-04-3] | 1000 | µg/mL | | |
| | Danitol | [39515-41-8] | 1000 | µg/mL | | |
| | Phenothrin | [26002-80-2] | 1000 | µg/mL | | |
| | Acrinathrin | [101007-06-1] | 1000 | µg/mL | | |
| | Cyhalothrin | [68085-85-8] | 1000 | µg/mL | | |
| | Permethrin | [52645-53-1] | 1000 | µg/mL | | |
| | beta-Cyfluthrin | [68359-37-5] | 1000 | µg/mL | | |
| | Cypermethrin | [52315-07-8] | 1000 | µg/mL | | |
| | Flucythrinate | [70124-77-5] | 1000 | µg/mL | | |
| | Fenvalerate | [51630-58-1] | 1000 | µg/mL | | |
| | Deltamethrin | [52918-63-5] | 1000 | µg/mL | | |
| | Tefluthrin | [79538-32-2] | 1000 | µg/mL | | |
| | Diphenylamine | [122-39-4] | 1000 | µg/mL | | |
| | Chlordimeform free base | [6164-98-3] | 1000 | µg/mL | | |
| | Trifluralin | [1582-09-8] | 1000 | µg/mL | | |
| | Dichlofluanid | [1085-98-9] | 1000 | µg/mL | | |
| | Fipronil | [120068-37-3] | 1000 | µg/mL | | |
| | Triadimefon | [43121-43-3] | 1000 | µg/mL | | |
| | Triadimenol | [55219-65-3] | 1000 | µg/mL | | |
| | Chlorfenapyr | [122453-73-0] | 1000 | µg/mL | | |
| | Flumetralin | [62924-70-3] | 1000 | µg/mL | | |
| | Dimepiperate | [61432-55-1] | 1000 | µg/mL | | |
| | Butralin | [33629-47-9] | 1000 | µg/mL | | |
| | Pendimethalin | [40487-42-1] | 1000 | µg/mL | | |
| | Quizalofop-ethyl | [76578-14-8] | 1000 | µg/mL | | |