

### **QuiKResponse - EDB/DBCP/TCP**

Art. ID	ERA-692QR
Unit	2 mL (ampoule)
Deliverydetails	No Dangerous Good

#### **Description**

One 2 mL flame-sealed ampoule yields in excess of 200 mL after dilution. Use with EPA method 8011. Each lot contains all analytes at .2 - 2.0 µg/L /// Analytes: 1,2-Dibromo-3-chloropropane (DBCP), 1,2-Dibromoethane (EDB), 1,2,3-Trichloropropane (TCP)