

USP Sterile Water System Suitability Set - 40 mL

Art. ID	ERA-18061
Unit	Set
Deliverydetails	No Dangerous Good /not restricted

Description

USP 643 Sterile Water System Suitability Set - Set of three vials includes (1) Reagent Water Rw, (1) 8.0 mg/L C from USP Sucrose Rs, (1) 8.0 mg/L C from USP 1,4-Benzoquinone Rss. Use for the following instruments: Sievers 800, 900, Analytik Jena, Shimadzu UV/Persulfate, & Teledyne Tekmar.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Reagent Water Rw	Reagent Water Rw					
USP Sucrose Rs	Sucrose	[57-50-1]	8	mg/L		
USP 1,4-Benzoquinone Rs	1,4-Benzochinone	[106-51-4]	8	mg/L		

s