

### **Ti-Alloy 6Al-2Mo-2Sn-4Zr (chips)**

Art. ID                      NIST-647  
Unit                         50 g  
Deliverydetails            No Dangerous Good /not restricted

#### Description

This Standard Reference Material (SRM®) is in the form of chips sized between 0.50 and 1.18 mm sieve openings (35 and 16 mesh). It is intended for use primarily in chemical methods analysis. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	5,88	%		
	Molybdenum (Mo)	[7439-98-7]	1,96	%		
	Tin (Sn)	[7440-31-5]	2,02	%		
	Zirconium (Zr)	[7440-67-7]	3,9	%		
	Iron (Fe)	[7439-89-6]	0,075	%		
	Carbon (C)	[7440-44-0]	0,006	%		