

**TT - 20.0 mg/L C from NIST KHP - 125 mL**

|                 |                   |
|-----------------|-------------------|
| Art. ID         | ERA-12903         |
| Unit            | 125 mL            |
| Deliverydetails | No Dangerous Good |

**Description**

TT- 20.0 mg/L C from NIST KHP -- 125 mL H<sub>3</sub>PO<sub>4</sub> - preserved NIST KHP standard appropriate for Teledyne Tekmar, OI or Shimadzu DO NOT USE FOR ANATEL 643 OR SIEVERS 400 6-month expiration