

## **Sievers Model 900 Calibration Kit**

|                 |                   |
|-----------------|-------------------|
| Art. ID         | ERA-19610         |
| Unit            | Kit               |
| Deliverydetails | No Dangerous Good |

### Description

Set includes eleven 40 mL vials, (1) blank, (1) each of 1.00, 5.00, 10.0, 25.0 and 50.0 mg/L C NIST KHP, (1) each of 1.00, 5.00, 10.0, 25.0 and 50.0 mg/L IC NIST NaHCO<sub>3</sub>. Replaces Sievers® Part STD 90000.