

Stainless Steel

Art. ID NCS HC28320
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Carbon (C) | [7440-44-0] | 0,013 | | | |
| | Silicon (Si) | [7440-21-3] | 0,464 | | | |
| | Manganese (Mn) | [7439-96-5] | 1,17 | | | |
| | Phosphorus (P) | [7723-14-0] | 0,033 | | | |
| | Sulfur (S) | [7704-34-9] | 0,0021 | | | |
| | Chromium (Cr) | [7440-47-3] | 19,67 | | | |
| | Nickel (Ni) | [7440-02-0] | 24,23 | | | |
| | Copper (Cu) | [7440-50-8] | 1,49 | | | |
| | Molybdenum (Mo) | [7439-98-7] | 4,47 | | | |
| | Vanadium (V) | [7440-62-2] | 0,086 | | | |
| | Titanium (Ti) | [7440-32-6] | 0,002 | | | |
| | Aluminium (Al) | [7429-90-5] | 0,024 | | | |
| | Cobalt (Co) | [7440-48-4] | 0,096 | | | |
| | Niobium (Nb) | [7440-03-1] | 0,042 | | | |
| | Tin (Sn) | [7440-31-5] | 0,0043 | | | |
| | Lead (Pb) | [7439-92-1] | 0,0001 | | | |
| | Arsenic (As) | [7440-38-2] | 0,0053 | | | |