

Refined Cast Iron

Art. ID	NIST-131h
Unit	100 g
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM) is intended primarily for use in validation of chemical and instrumental methods of analysis. It can be used to validate value assignment of in-house reference materials and, if necessary, to calibrate carbon/sulfur analyzers. A unit of NIST131h consists of a bottle containing 100 g of chips. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Sulfur (S)	[7704-34-9]	7,4	mg/kg		
Reference value	Carbon (C)	[7440-44-0]	7,8	mg/kg		