

## **Iodized Table Salt (Iodide)**

Art. ID	NIST-3530
Unit	1 bottle x 200 g
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is intended primarily for validation of methods for the determination of iodide in an iodized table salt or similar material. This SRM can also be used for quality assurance when assigning values to in-house reference materials. A unit of NIST-3530 consists of one bottle containing approximately 200 g of a commercially available iodized table salt sealed inside an aluminized polyethylene terephthalate (PET) pouch containing an oxygen absorber pack. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Iodine (as Iodide)		52,2 ± 4,2	mg/kg		