

Custom Solution - 4 components in 2% HNO₃

Art. ID LM24-CUS-70081
Unit 2 x 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	150	µg/mL		
	Cadmium (Cd)	[7440-43-9]	50	µg/mL		
	Thallium (Tl)	[7440-28-0]	80	µg/mL		
	Lead (Pb)	[7439-92-1]	50	µg/mL		