

### **Custom Solution - 5 components in 5% HNO<sub>3</sub>**

Art. ID LM24-CUS-70082  
Unit 2 x 25 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 II E2

#### **Description**

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	140	mg/L		
	Lithium (Li)	[7439-93-2]	55	mg/L		
	Molybdenum (Mo)	[7439-98-7]	300	mg/L		
	Chromium (Cr)	[7440-47-3]	1100	mg/L		
	Copper (Cu)	[7440-50-8]	300	mg/L		