

Custom Solution - 9 components in 5% HCl

Art. ID LM24-CUS-70083
Unit 2 x 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	120	mg/L		
	Iridium (Ir)	[7439-88-5]	10	mg/L		
	Gold (Au)	[7440-57-5]	10	mg/L		
	Tin (Sn)	[7440-31-5]	600	mg/L		
	Osmium (Os)	[7440-04-2]	10	mg/L		
	Palladium (Pd)	[7440-05-3]	10	mg/L		
	Platinum (Pt)	[7440-06-4]	10	mg/L		
	Rhodium (Rh)	[7440-16-6]	10	mg/L		
	Ruthenium (Ru)	[7440-18-8]	10	mg/L		