

Custom Solution - 5 components in 5% HNO₃

Art. ID LM24-CUS-70084
Unit 2 x 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	150	mg/L		
	Selenium (Se)	[7782-49-2]	150	mg/L		
	Cobalt (Co)	[7440-48-4]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	200	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		