

Trace Elements in Water (lot 1219) - TMDA-64.3 sold out

Art. ID NW-TMDA-64.3
Unit 250 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	293 ± 25	µg/L		
	Antimony (Sb)	[7440-36-0]	124 ± 10	µg/L		
	Arsenic (As)	[7440-38-2]	163 ± 13	µg/L		
	Barium (Ba)	[7440-39-3]	287 ± 18	µg/L		
	Beryllium (Be)	[7440-41-7]	157 ± 13	µg/L		
	Bismuth (Bi)	[7440-69-9]	117 ± 24	µg/L		
	Boron (B)	[7440-42-8]	277 ± 30	µg/L		
	Cadmium (Cd)	[7440-43-9]	258 ± 20	µg/L		
	Chromium (Cr)	[7440-47-3]	283 ± 18	µg/L		
	Cobalt (Co)	[7440-48-4]	248 ± 16	µg/L		
	Copper (Cu)	[7440-50-8]	261 ± 18	µg/L		
	Gallium (Ga)	[7440-55-3]	50,3 ± 3,4	µg/L		
	Iron (Fe)	[7439-89-6]	299 ± 20	µg/L		
	Lead	[7439-92-1]	279 ± 22	µg/L		
	Lithium (Li)	[7439-93-2]	144 ± 14	µg/L		
	Manganese (Mn)	[7439-96-5]	292 ± 16	µg/L		
	Molybdenum (Mo)	[7439-98-7]	287 ± 20	µg/L		
	Nickel (Ni)	[7440-02-0]	251 ± 17	µg/L		
	Rubidium (Rb)	[7440-17-7]	30,8 ± 2,3	µg/L		
	Selenium (Se)	[7782-49-2]	154 ± 19	µg/L		
	Silver (Ag)	[7440-22-4]	12,4 ± 1,4	µg/L		
	Strontium (Sr)	[7440-24-6]	630 ± 33	µg/L		
	Thallium (Tl)	[7440-28-0]	140 ± 15	µg/L		
	Tin	[7440-31-5]	281 ± 29	µg/L		
	Titanium (Ti)	[7440-32-6]	124 ± 7	µg/L		
	Tungsten (W)	[7440-33-7]	0,382 ± 0,068	µg/L		
	Uranium (U)	[7440-61-1]	134 ± 11	µg/L		
	Vanadium (V)	[7440-62-2]	280 ± 18	µg/L		
	Zinc (Zn)	[7440-66-6]	320 ± 26	µg/L		