

ARCONIC - Aluminium base , Type : Cu, Class: RM, cast, Disc 64 mm x 25 mm

Art. ID ALCOA-CU-2
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 2,5i | % | | |
| | Iron (Fe) | [7439-89-6] | 0,25i | % | | |
| | Manganese (Mn) | [7439-96-5] | <0,0050 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,15i | % | | |