

**ARCONIC - Aluminium base , Type : Cu, Class: RM, cast, Disc 64 mm x 25 mm**

Art. ID ALCOA-CU-2  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	2,5i	%		
	Iron (Fe)	[7439-89-6]	0,25i	%		
	Manganese (Mn)	[7439-96-5]	<0,0050	%		
	Silicon (Si)	[7440-21-3]	0,15i	%		