

ARCONIC - Aluminium base , Type : Si, Class: RM, cast, Disc 64 mm x 25 mm

Art. ID ALCOA-SI-3
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Iron (Fe) | [7439-89-6] | ~0,5 | % | | |
| | Silicon (Si) | [7440-21-3] | 1,6i | % | | |