

**ARCONIC - Aluminium base 7047, Class: RM, Alloy: 7047, cast, typical analysis, ppm values not shown, Disc 64 mm x 25 mm**

Art. ID ALCOA-SS-7047  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,368	%		
	Chromium (Cr)	[7440-47-3]	0,011	%		
	Copper (Cu)	[7440-50-8]	0,024	%		
	Iron (Fe)	[7439-89-6]	0,051	%		
	Gallium (Ga)	[7440-55-3]	0,01	%		
	Magnesium (Mg)	[7439-95-4]	1,54	%		
	Manganese (Mn)	[7439-96-5]	0,021	%		
	Silicon (Si)	[7440-21-3]	0,042	%		
	Titanium (Ti)	[7440-32-6]	0,031	%		
	Vanadium (V)	[7440-62-2]	0,011	%		
	Zinc (Zn)	[7440-66-6]	7,48	%		
	Zirconium (Zr)	[7440-67-7]	0,106	%		