

ARCONIC - Aluminium base 8011, Class: RM, Alloy: 8011, cast, typical analysis, ppm values not shown, Disc 64 mm x 25 mm

Art. ID ALCOA-SS-8011
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 0,0218 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,699 | % | | |
| | Gallium (Ga) | [7440-55-3] | 0,0147 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,0217 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,536 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,0394 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,0129 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,0201 | % | | |