

**ARCONIC - Aluminium base 8011, Class: RM, Alloy: 8011, cast, typical analysis, ppm values not shown, Disc 64 mm x 25 mm**

Art. ID ALCOA-SS-8011  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	0,0218	%		
	Iron (Fe)	[7439-89-6]	0,699	%		
	Gallium (Ga)	[7440-55-3]	0,0147	%		
	Manganese (Mn)	[7439-96-5]	0,0217	%		
	Silicon (Si)	[7440-21-3]	0,536	%		
	Titanium (Ti)	[7440-32-6]	0,0394	%		
	Vanadium (V)	[7440-62-2]	0,0129	%		
	Zinc (Zn)	[7440-66-6]	0,0201	%		